

---

# **APS\_BlueSky\_tools Documentation**

*Release 0.0.40+28.g0822f71.dirty*

**Pete R. Jemian**

**Jan 03, 2019**



---

## Contents:

---

|          |                                |           |
|----------|--------------------------------|-----------|
| <b>1</b> | <b>Package Information</b>     | <b>3</b>  |
| 1.1      | Applications . . . . .         | 3         |
| 1.2      | bluesky_snapshot . . . . .     | 3         |
| 1.3      | Examples . . . . .             | 7         |
| 1.4      | Callbacks . . . . .            | 9         |
| 1.5      | Devices . . . . .              | 10        |
| 1.6      | File Writers . . . . .         | 23        |
| 1.7      | Plans . . . . .                | 26        |
| 1.8      | Signals . . . . .              | 30        |
| 1.9      | Suspenders . . . . .           | 31        |
| 1.10     | Utilities . . . . .            | 31        |
| 1.11     | synApps busy record . . . . .  | 34        |
| 1.12     | synApps sscan record . . . . . | 34        |
| 1.13     | synApps swait record . . . . . | 35        |
| <b>2</b> | <b>Indices and tables</b>      | <b>37</b> |
|          | <b>Python Module Index</b>     | <b>39</b> |



Various Python tools for use with BlueSky at the APS

- <http://nsls-ii.github.io/>
- <https://github.com/NSLS-II/bluesky>



# CHAPTER 1

---

## Package Information

---

**author**

Pete R. Jemian

**email** [jemian@anl.gov](mailto:jemian@anl.gov)

**copyright** 2017-2019, Pete R. Jemian

**license** ANL OPEN SOURCE LICENSE (see LICENSE file)

**documentation** [https://APS\\_BlueSky\\_tools.readthedocs.io](https://APS_BlueSky_tools.readthedocs.io)

**source** [https://github.com/BCDA-APS/APS\\_BlueSky\\_tools](https://github.com/BCDA-APS/APS_BlueSky_tools)

## 1.1 Applications

There are two command-line applications provided by APS\_BlueSky\_tools:

| application purpose                            |   |
|--|---|
| <a href="#">aps_blueSky_tools_plan_catalog</a> | summary list of all scans in the databroker                             |
| <a href="#">bluesky_snapshot</a>               | Take a snapshot of a list of EPICS PVs and record it in the databroker. |

## 1.2 bluesky\_snapshot

Take a snapshot of a list of EPICS PVs and record it in the databroker. Retrieve (and display) that snapshot later using `APS_BlueSky_tools.callbacks.SnapshotReport`.

### 1.2.1 Example - command line

Before using the command-line interface, find out what the `bluesky_snapshot` expects:

```
$ bluesky_snapshot -h
usage: bluesky_snapshot [-h] [-b BROKER_CONFIG] [-m METADATA_SPEC] [-r] [-v]
                        EPICS_PV [EPICS_PV ...]

record a snapshot of some PVs using Bluesky, ophyd, and databroker
version=0.0.40+26.g323cd35

positional arguments:
  EPICS_PV            EPICS PV name

optional arguments:
  -h, --help          show this help message and exit
  -b BROKER_CONFIG   YAML configuration for databroker, default:
                     mongodb_config
  -m METADATA_SPEC   --metadata METADATA_SPEC
                     additional metadata, enclose in quotes, such as -m
                     "purpose=just tuned, situation=routine"
  -r, --report        suppress snapshot report
  -v, --version       show program's version number and exit
```

The help does not tell you that the default for BROKER\_CONFIG is “mongodb\_config”, a YAML file in one of the default locations where the databroker expects to find it. That’s what we have.

We want to snapshot just a couple PVs to show basic use. Here are their current values:

```
$ caget prj:IOC_CPU_LOAD prj:SYS_CPU_LOAD
prj:IOC_CPU_LOAD          0.900851
prj:SYS_CPU_LOAD          4.50426
```

Here’s the snapshot (we’ll also set a metadata that says this is an example):

```
$ bluesky_snapshot prj:IOC_CPU_LOAD prj:SYS_CPU_LOAD -m "purpose=example"
=====
snapshot: 2019-01-03 17:02:42.922197
=====

hints: {}
hostname: mint-vm
iso8601: 2019-01-03 17:02:42.922197
login_id: mintadmin@mint-vm
plan_description: archive snapshot of ophyd Signals (usually EPICS PVs)
plan_name: snapshot
plan_type: generator
purpose: example
scan_id: 1
software_versions: {'python': '3.6.6 |Anaconda custom (64-bit)| (default, Jun 28 2018,
                   ↵ 17:14:51) \n[GCC 7.2.0]', 'PyEpics': '3.3.1', 'bluesky': '1.4.1', 'ophyd': '1.3.0',
                   ↵ 'databroker': '0.11.3', 'APS_Bluesky_Tools': '0.0.40+26.g323cd35.dirty'}
time: 1546556562.9231327
uid: 98a86a91-d41e-4965-a048-afa5b982a17c
username: mintadmin

===== ===== ===== =====
timestamp           source name      value
===== ===== ===== =====
2019-01-03 17:02:33.930067 PV    prj:IOC_CPU_LOAD 0.8007421685989062
2019-01-03 17:02:33.930069 PV    prj:SYS_CPU_LOAD 10.309472772459404
```

(continues on next page)

(continued from previous page)

```
=====
exit_status: success
num_events: {'primary': 1}
run_start: 98a86a91-d41e-4965-a048-afa5b982a17c
time: 1546556563.1087885
uid: 026fa69c-45b7-4b45-a3b3-266aadbf7176
```

We have a second IOC (`gov`) that has the same PVs. Let's get them, too.:

```
$ bluesky_snapshot {gov,otz}:{IOC,SYS}_CPU_LOAD -m "purpose=this is an example,_
˓→example=example 2"

=====
snapshot: 2018-12-20 18:21:53.371995
=====

example: example 2
hints: {}
iso8601: 2018-12-20 18:21:53.371995
plan_description: archive snapshot of ophyd Signals (usually EPICS PVs)
plan_name: snapshot
plan_type: generator
purpose: this is an example
scan_id: 1
software_versions: {'python': '3.6.2 |Continuum Analytics, Inc.| (default, Jul 20_
˓→2017, 13:51:32) \n[GCC 4.4.7 20120313 (Red Hat 4.4.7-1)]', 'PyEpics': '3.3.1',
˓→'bluesky': '1.4.1', 'ophyd': '1.3.0', 'databroker': '0.11.3', 'APS_Bluesky_Tools':
˓→'0.0.37'}
time: 1545351713.3727024
uid: d5e15ba3-0393-4df3-8217-1b72d82b5cf9

=====
timestamp          source name           value
=====
2018-12-20 18:21:45.488033 PV      gov:IOC_CPU_LOAD 0.22522293126578166
2018-12-20 18:21:45.488035 PV      gov:SYS_CPU_LOAD 10.335244804189122
2018-12-20 18:21:46.910976 PV      otz:IOC_CPU_LOAD 0.10009633509509736
2018-12-20 18:21:46.910973 PV      otz:SYS_CPU_LOAD 11.360899731293234
=====

exit_status: success
num_events: {'primary': 1}
run_start: d5e15ba3-0393-4df3-8217-1b72d82b5cf9
time: 1545351713.3957422
uid: e033cd99-dcac-4b56-848c-62eede1e4d77
```

You can log text and arrays, too.:

```
$ bluesky_snapshot {gov,otz}:{iso8601,HOSTNAME,{IOC,SYS}_CPU_LOAD} compress \
-m "purpose=this is an example, example=example 2, look=can snapshot text and_
˓→arrays too, note=no commas in metadata"

=====
snapshot: 2018-12-20 18:28:28.825551
=====
```

(continues on next page)

(continued from previous page)

```

example: example 2
hints: {}
iso8601: 2018-12-20 18:28:28.825551
look: can snapshot text and arrays too
note: no commas in metadata
plan_description: archive snapshot of ophyd Signals (usually EPICS PVs)
plan_name: snapshot
plan_type: generator
purpose: this is an example
scan_id: 1
software_versions: {'python': '3.6.2 |Continuum Analytics, Inc.| (default, Jul 20_
˓→2017, 13:51:32) \n[GCC 4.4.7 20120313 (Red Hat 4.4.7-1)]', 'PyEpics': '3.3.1',
˓→'bluesky': '1.4.1', 'ophyd': '1.3.0', 'databroker': '0.11.3', 'APS_Bluessky_Tools':
˓→'0.0.37'}
time: 1545352108.8262713
uid: 7e77708e-9169-45ab-b2b6-4e31534d980a

=====
timestamp          source name        value
=====
2018-12-20 18:24:34.220028 PV    compress      [0.1, 0.2, 0.3]
2018-12-13 14:49:53.121188 PV    gov:HOSTNAME otz.aps.anl.gov
2018-12-20 18:28:25.093941 PV    gov:IOC_CPU_LOAD 0.1501490058473918
2018-12-20 18:28:25.093943 PV    gov:SYS_CPU_LOAD 10.360270546421546
2018-12-20 18:28:28.817630 PV    gov:iso8601   2018-12-20T18:28:28
2018-12-13 14:49:53.135016 PV    otz:HOSTNAME otz.aps.anl.gov
2018-12-20 18:28:26.525208 PV    otz:IOC_CPU_LOAD 0.10009727705620367
2018-12-20 18:28:26.525190 PV    otz:SYS_CPU_LOAD 12.937574161543873
2018-12-20 18:28:28.830285 PV    otz:iso8601   2018-12-20T18:28:28
=====

exit_status: success
num_events: {'primary': 1}
run_start: 7e77708e-9169-45ab-b2b6-4e31534d980a
time: 1545352108.8656788
uid: 0de0ec62-504e-4dbc-ad08-2507d4ed44f9

```

## 1.2.2 Source code documentation

record a snapshot of some PVs using Bluesky, ophyd, and databroker

**USAGE:**

```

(base) user@hostname .../pwd $ bluesky_snapshot -h
usage: bluesky_snapshot [-h] [-b BROKER_CONFIG] [-m METADATA_SPEC] [-r] [-v]
                        EPICS_PV [EPICS_PV ...]

record a snapshot of some PVs using Bluesky, ophyd, and databroker
version=0.0.40+26.g323cd35

positional arguments:
  EPICS_PV            EPICS PV name

optional arguments:
  -h, --help           show this help message and exit

```

(continues on next page)

(continued from previous page)

|  |   |
|--|---|
| -b BROKER_CONFIG                           | YAML configuration for databroker, default:<br>mongodb_config                                 |
| -m METADATA_SPEC, --metadata METADATA_SPEC | additional metadata, enclose in quotes, such as -m<br>"purpose=just tuned, situation=routine" |
| -r, --report                               | suppress snapshot report  |
| -v, --version                              | show program's version number and exit  |

`APS_BlueSky_tools.snapshot.get_args()`  
get command line arguments

`APS_BlueSky_tools.snapshot.snapshot_cli()`  
given a list of PVs on the command line, snapshot and print report

EXAMPLES:

```
snapshot.py pv1 [more pvs ...]
snapshot.py `cat pvlist.txt`
```

Note that these are equivalent:

```
snapshot.py rpi5bf5:0:humidity rpi5bf5:0:temperature
snapshot.py rpi5bf5:0:{humidity,temperature}
```

## 1.3 Examples

- *Example: plan\_catalog()*
- *Example: specfile\_example()*
- *Example: nscan()*
- *Example: TuneAxis()*
- *Source Code Documentation*
- *Downloads*

### 1.3.1 Example: `plan_catalog()`

The APS\_BlueSky\_tools package provides an executable that can be used to display a summary of all the scans in the database. The executable wraps the demo function: `plan_catalog()`. It is for demonstration purposes only (since it does not filter the output to any specific subset of scans).

The output is a table, formatted as restructured text, with these columns:

**date/time** The date and time the scan was started.

**short\_uid** The first characters of the scan's UUID (unique identifier).

**id** The scan number. (User has control of this and could reset the counter for the next scan.)

**plan** Name of the plan that initiated this scan.

**args** Arguments to the plan that initiated this scan.

This is run as a linux console command:

```
aps_blueSky_tools_plan_catalog | tee out.txt
```

The full output is almost a thousand lines. Here are the first few lines:

```
1 ====== ===== ===== ====== 
2 date/time      short_uid id   plan          args
3 
4 2017-10-26 11:21:28 3fe59011 1   scan          detectors=['noisy'],
5   num=219, motor=['m1'], start=-1.5, stop=-0.5, per_step=None
6 
7 2017-10-26 11:21:42 25b4c903 2   scan          detectors=['noisy'],
8   num=219, motor=['m1'], start=-1.5, stop=-0.5, per_step=None
9 
10 2017-10-26 11:22:08 3953e8e0 3   scan         detectors=['noisy'],
11   num=219, motor=['m1'], start=-1.5, stop=-0.5, per_step=None
12 
13 2017-10-26 11:22:22 f24bf2cc 4   scan         detectors=['noisy'],
14   num=219, motor=['m1'], start=-1.5, stop=-0.5, per_step=None
15 
16 2017-10-26 11:22:37 44b751d2 5   scan         detectors=['noisy'],
17   num=219, motor=['m1'], start=-1.5, stop=-0.5, per_step=None
18 
19 2017-10-26 11:22:50 4e3741f5 6   scan         detectors=['noisy'],
20   num=219, motor=['m1'], start=-1.5, stop=-0.5, per_step=None
21 
22 2017-10-26 11:24:33 a83df5d4 7   scan         detectors=['synthetic_',
23   pseudovoigt'], num=219, motor=['m1'], start=-2, stop=0, per_step=None
24
```

### 1.3.2 Example: `specfile_example()`

We'll use a Jupyter notebook to demonstrate the `specfile_example()` that writes one or more scans to a SPEC data file. Follow here: [https://github.com/BCDA-APS/APS\\_BlueSky\\_tools/blob/master/docs/source/resources/demo\\_specfile\\_example.ipynb](https://github.com/BCDA-APS/APS_BlueSky_tools/blob/master/docs/source/resources/demo_specfile_example.ipynb)

### 1.3.3 Example: `nscan()`

We'll use a Jupyter notebook to demonstrate the `nscan()` plan. An *nscan* is used to scan two or more axes together, such as a  $\theta$ - $2\theta$  diffractometer scan. Follow here: [https://github.com/BCDA-APS/APS\\_BlueSky\\_tools/blob/master/docs/source/resources/demo\\_nscan.ipynb](https://github.com/BCDA-APS/APS_BlueSky_tools/blob/master/docs/source/resources/demo_nscan.ipynb)

### 1.3.4 Example: `TuneAxis()`

We'll use a Jupyter notebook to demonstrate the `TuneAxis()` support that provides custom alignment of a signal against an axis. Follow here: [https://github.com/BCDA-APS/APS\\_BlueSky\\_tools/blob/master/docs/source/resources/demo\\_tuneaxis.ipynb](https://github.com/BCDA-APS/APS_BlueSky_tools/blob/master/docs/source/resources/demo_tuneaxis.ipynb)

### 1.3.5 Source Code Documentation

demonstrate BlueSky callbacks

---

|  |   |
|--|---|
| <code>plan_catalog(db)</code>                      | make a table of all scans known in the databroker     |
| <code>specfile_example(headers[, filename])</code> | write one or more headers (scans) to a SPEC data file |

---

`APS_BlueSky_tools.examples.main()`  
summary list of all scans in the databroker

`aps_blueSky_tools_plan_catalog` command-line application

This can be unwieldy if there are many scans in the databroker. Consider it as a demo program rather than for general, long-term use.

`APS_BlueSky_tools.examples.plan_catalog(db)`  
make a table of all scans known in the databroker

Example:

```
from APS_BlueSky_tools.examples import plan_catalog
plan_catalog(db)
```

`APS_BlueSky_tools.examples.specfile_example(headers, filename='test_specdata.txt')`  
write one or more headers (scans) to a SPEC data file

### 1.3.6 Downloads

The jupyter notebook and files related to this section may be downloaded from the following table.

- `plan_catalog.txt`
- jupyter notebook: `demo_nscan`
- jupyter notebook: `demo_tuneaxis`
- jupyter notebook: `demo_specfile_example`
  - `spec1.dat`
  - `spec2.dat`
  - `spec3.dat`
  - `spec_tunes.dat`
  - `test_specdata.txt`

## 1.4 Callbacks

Callbacks that might be useful at the APS using BlueSky

---

|   |  |
|---|--|
| <code>document_contents_callback(key, doc)</code> | prints document contents – use for diagnosing a document stream        |
| <code>DocumentCollectorCallback()</code>          | BlueSky callback to collect <i>all</i> documents from most-recent plan |

---

Continued on next page

Table 2 – continued from previous page

|  |  |
|--|--|
| <code>SnapshotReport(*args, **kwargs)</code> | show the data from a <code>APS_BlueSky_Tools.plans.snapshot()</code> |
|--|--|

## FILE WRITER CALLBACK

see `SpecWriterCallback()`

**class** `APS_BlueSky_tools.callbacks.DocumentCollectorCallback`  
BlueSky callback to collect *all* documents from most-recent plan

Will reset when it receives a *start* document.

EXAMPLE:

```
from APS_BlueSky_tools.callbacks import DocumentCollector
doc_collector = DocumentCollectorCallback()
RE.subscribe(doc_collector.receiver)
...
RE(some_plan())
print(doc_collector.uids)
print(doc_collector.documents["stop"])
```

**receiver** (*key, document*)

keep all documents from recent plan in memory

**class** `APS_BlueSky_tools.callbacks.SnapshotReport (*args, **kwargs)`  
show the data from a `APS_BlueSky_Tools.plans.snapshot()`

Find most recent snapshot between certain dates:

```
headers = db(plan_name="snapshot", since="2018-12-15", until="2018-12-21")
h = list(headers)[0]           # pick the first one, it's the most recent
APS_BlueSky_Tools.callbacks.SnapshotReport().print_report(h)
```

Use as callback to a snapshot plan:

```
RE(
    APS_BlueSky_Tools.plans.snapshot(ophyd_objects_list),
    APS_BlueSky_Tools.callbacks.SnapshotReport()
)
```

**descriptor** (*doc*)

**special case:** the data is both in the descriptor AND the event docs due to the way our plan created it

**print\_report** (*header*)

simplify the job of writing our custom data table

method: play the entire document stream through this callback

`APS_BlueSky_tools.callbacks.document_contents_callback(key, doc)`  
prints document contents – use for diagnosing a document stream

## 1.5 Devices

(ophyd) Devices that might be useful at the APS using BlueSky

APS GENERAL SUPPORT

---

|  |  |
|--|--|
| <code>ApsMachineParametersDevice(*args, **kwargs)</code>       | common operational parameters of the APS of general interest |
| <code>ApsPssShutter(*args, **kwargs)</code>                    | APS PSS shutter  |
| <code>ApsPssShutterWithStatus(prefix, state_pv, ...)</code>    | APS PSS shutter with separate status PV                      |
| <code>SimulatedApsPssShutterWithStatus(*args, **kwargs)</code> | Simulated APS PSS shutter                                    |

---

## AREA DETECTOR SUPPORT

---

|  |  |
|--|--|
| <code>AD_setup_FrameType(prefix[, scheme])</code>  | configure so frames are identified & handled by type (dark, white, or image) |
| <code>AD_warmed_up(detector)</code>                | Has area detector pushed an NDarray to the HDF5 plugin? True or False        |
| <code>AD_EpicsHdf5FileName(*args, **kwargs)</code> | custom class to define image file name from EPICS                            |

---

## DETECTOR / SCALER SUPPORT

---

|  |   |
|--|---|
| <code>use_EPICS_scaler_channels(scaler)</code> | configure scaler for only the channels with names assigned in EPICS |
|--|---|

---

## MOTORS, POSITIONERS, AXES, ...

---

|   |   |
|---|---|
| <code>AxisTunerException</code>                     | Exception during execution of <code>AxisTunerBase</code> subclass                                   |
| <code>AxisTunerMixin(*args, **kwargs)</code>        | Mixin class to provide tuning capabilities for an axis  |
| <code>EpicsDescriptionMixin(*args, **kwargs)</code> | add a record's description field to a Device, such as <code>EpicsMotor</code>                       |
| <code>EpicsMotorDialMixin(*args, **kwargs)</code>   | add motor record's dial coordinate fields to Device   |
| <code>EpicsMotorLimitsMixin(*args, **kwargs)</code> | add motor record HLM & LLM fields & compatibility <code>get_lim()</code> and <code>set_lim()</code> |
| <code>EpicsMotorRawMixin(*args, **kwargs)</code>    | add motor record's raw coordinate fields to Device  |
| <code>EpicsMotorServoMixin(*args, **kwargs)</code>  | add motor record's servo loop controls to Device  |
| <code>EpicsMotorShutter(*args, **kwargs)</code>     | a shutter, implemented with an EPICS motor moved between two positions                              |
| <code>EpicsOnOffShutter(*args, **kwargs)</code>     | a shutter, implemented with an EPICS PV moved between two positions                                 |

---

## SHUTTERS

---

|   |  |
|---|--|
| <code>ApsPssShutter(*args, **kwargs)</code>                 | APS PSS shutter  |
| <code>ApsPssShutterWithStatus(prefix, state_pv, ...)</code> | APS PSS shutter with separate status PV                                |
| <code>EpicsMotorShutter(*args, **kwargs)</code>             | a shutter, implemented with an EPICS motor moved between two positions |
| <code>EpicsOnOffShutter(*args, **kwargs)</code>             | a shutter, implemented with an EPICS PV moved between two positions    |

---

## synApps records

---

|   |   |
|---|---|
| <code>busyRecord(*args, **kwargs)</code>  |   |
| <code>sscanRecord(*args, **kwargs)</code> | EPICS synApps sscan record: used as \$(P):scan(N) |

---

Continued on next page

Table 8 – continued from previous page

|  |   |
|--|---|
| <code>sscanDevice(*args, **kwargs)</code>                        | synApps XXX IOC setup of sscan records: \$(P):scan\$(N) |
| <code>sawaitRecord(*args, **kwargs)</code>                       | synApps sawait record: used as \$(P):userCalc\$(N)      |
| <code>sawait_setup_random_number(sawait, **kw)</code>            | setup sawait record to generate random numbers          |
| <code>sawait_setup_gaussian(sawait, motor[, center, ...])</code> | setup sawait for noisy Gaussian                         |
| <code>sawait_setup_lorentzian(sawait, motor[, ...])</code>       | setup sawait record for noisy Lorentzian                |
| <code>sawait_setup_incremter(sawait[, scan, limit])</code>       | setup sawait record as an incrementer                   |
| <code>userCalcsDevice(*args, **kwargs)</code>                    | synApps XXX IOC setup of userCalcs: \$(P):userCalc\$(N) |

## OTHER SUPPORT

|   |   |
|---|---|
| <code>DualPF4FilterBox(*args, **kwargs)</code>      | Dual Xia PF4 filter boxes using support from synApps (using Al, Ti foils) |
| <code>EpicsDescriptionMixin(*args, **kwargs)</code> | add a record's description field to a Device, such as EpicsMotor          |
| <code>ProcedureRegistry(*args, **kwargs)</code>     | Procedure Registry: run a blocking function in a thread                   |

## Internal routines

|   |  |
|---|--|
| <code>ApsOperatorMessagesDevice(*args, **kwargs)</code> | general messages from the APS main control room        |
| <code>DeviceMixinBase(*args, **kwargs)</code>           | Base class for APS_Bluessky_tools Device mixin classes |

`class APS_BlueSky_tools.devices.AD_EpicsHdf5FileName (*args, **kwargs)`  
custom class to define image file name from EPICS

**Caution:** *Caveat emptor* applies here. You assume expertise!

Replace standard Bluesky algorithm where file names are defined as UUID strings, virtually guaranteeing that no existing images files will ever be overwritten.

Also, this method decouples the data files from the databroker, which needs the files to be named by UUID.

|  |   |
|--|---|
| <code>make_filename()</code>                             | overrides default behavior: Get info from EPICS HDF5 plugin.  |
| <code>generate_datum(key, timestamp, datum_kwarg)</code> | Generate a uid and cache it with its key for later insertion. |
| <code>get_frames_per_point()</code>                      | overrides default behavior                                    |
| <code>stage()</code>                                     | overrides default behavior                                    |

To allow users to control the file **name**, we override the `make_filename()` method here and we need to override some intervening classes.

To allow users to control the file **number**, we override the `stage()` method here and triple-comment out that line, and bring in sections from the methods we are replacing here.

The image file name is set in `FileStoreBase.make_filename()` from `ophyd.areadetector.filestore_mixins`. This is called (during device staging) from `FileStoreBase.stage()`

EXAMPLE:

To use this custom class, we need to connect it to some intervening structure. Here are the steps:

1. override default file naming
2. use to make your custom iterative writer
3. use to make your custom HDF5 plugin
4. use to make your custom AD support

imports:

```
from bluesky import RunEngine, plans as bp
from ophyd.areadetector import SimDetector, SingleTrigger
from ophyd.areadetector import ADComponent, ImagePlugin, SimDetectorCam
from ophyd.areadetector import HDF5Plugin
from ophyd.areadetector.filestore_mixins import FileStoreIterativeWrite
```

override default file naming:

```
from APS_BlueSky_tools.devices import AD_EpicsHdf5FileName
```

make a custom iterative writer:

```
class myHdf5EpicsIterativeWriter(AD_EpicsHdf5FileName, FileStoreIterativeWrite):
    pass
```

make a custom HDF5 plugin:

```
class myHDF5FileNames(HDF5Plugin, myHdf5EpicsIterativeWriter): pass
```

define support for the detector (simulated detector here):

```
class MySimDetector(SingleTrigger, SimDetector):
    '''SimDetector with HDF5 file names specified by EPICS'''

    cam = ADComponent(SimDetectorCam, "cam1:")
    image = ADComponent(ImagePlugin, "image1:")

    hdf1 = ADComponent(
        myHDF5FileNames,
        suffix = "HDF1:",
        root = "/",
        write_path_template = "/",
    )
```

create an instance of the detector:

```
simdet = MySimDetector("13SIM1:", name="simdet")
if hasattr(simdet.hdf1.stage_sigs, "array_counter"):
    # remove this so array counter is not set to zero each staging
    del simdet.hdf1.stage_sigs["array_counter"]
simdet.hdf1.stage_sigs["file_template"] = '%s%s_%3.3d.h5'
```

setup the file names using the EPICS HDF5 plugin:

```
simdet.hdf1.file_path.put("/tmp/simdet_demo/")      # ! ALWAYS end with a "/" !
simdet.hdf1.file_name.put("test")
simdet.hdf1.array_counter.put(0)
```

If you have not already, create a bluesky RunEngine:

```
RE = RunEngine({})
```

take an image:

```
RE(bp.count([simdet]))
```

## INTERNAL METHODS

**generate\_datum**(*key, timestamp, datum\_kwargs*)

Generate a uid and cache it with its key for later insertion.

**get\_frames\_per\_point**()

overrides default behavior

**make\_filename**()

overrides default behavior: Get info from EPICS HDF5 plugin.

**stage**()

overrides default behavior

Set EPICS items before device is staged, then copy EPICS naming template (and other items) to ophyd after staging.

`APS_BlueSky_tools.devices.AD_setup_FrameType(prefix, scheme='NeXus')`  
configure so frames are identified & handled by type (dark, white, or image)

## PARAMETERS

*prefix* (str) : EPICS PV prefix of area detector, such as “13SIM1:” *scheme* (str) : any key in the *AD\_FrameType\_schemes* dictionary

This routine prepares the EPICS Area Detector to identify frames by image type for handling by clients, such as the HDF5 file writing plugin. With the HDF5 plugin, the *FrameType* PV is added to the NDAttributes and then used in the layout file to direct the acquired frame to the chosen dataset. The *FrameType* PV value provides the HDF5 address to be used.

To use a different scheme than the defaults, add a new key to the *AD\_FrameType\_schemes* dictionary, defining storage values for the fields of the EPICS *mbbo* record that you will be using.

see: [https://github.com/BCDA-APS/use\\_bluesky/blob/master/notebooks/images\\_darks\\_flats.ipynb](https://github.com/BCDA-APS/use_bluesky/blob/master/notebooks/images_darks_flats.ipynb)

## EXAMPLE:

```
AD_setup_FrameType("2bmbPG3:", scheme="DataExchange")
```

- Call this function *before* creating the ophyd area detector object
- use lower-level PyEpics interface

`APS_BlueSky_tools.devices.AD_warmed_up(detector)`

Has area detector pushed an NDarray to the HDF5 plugin? True or False

Works around an observed issue: #598 <https://github.com/NSLS-II/ophyd/issues/598#issuecomment-414311372>

If detector IOC has just been started and has not yet taken an image with the HDF5 plugin, then a `TimeoutError` will occur as the HDF5 plugin “Capture” is set to 1 (Start). In such case, first acquire at least one image with the HDF5 plugin enabled.

---

**class** APS\_BlueSky\_tools.devices.**ApsBssUserInfoDevice**(\*args, \*\*kwargs)  
provide current experiment info from the APS BSS

BSS: Beamtime Scheduling System

EXAMPLE:

```
bss_user_info = ApsBssUserInfoDevice(  
    "9id_bss:",  
    name="bss_user_info")  
sd.baseline.append(bss_user_info)
```

**class** APS\_BlueSky\_tools.devices.**ApsMachineParametersDevice**(\*args, \*\*kwargs)  
common operational parameters of the APS of general interest

EXAMPLE:

```
import APS_BlueSky_tools.devices as APS_devices  
APS = APS_devices.ApsMachineParametersDevice(name="APS")  
aps_current = APS.current  
  
# make sure these values are logged at start and stop of every scan  
sd.baseline.append(APS)  
# record storage ring current as secondary stream during scans  
# name: aps_current_monitor  
# db[-1].table("aps_current_monitor")  
sd.monitors.append(aps_current)
```

The *sd.baseline* and *sd.monitors* usage relies on this global setup:

```
from bluesky import SupplementalData sd = SupplementalData() RE.preprocessors.append(sd)
```

---

|                         |  |
|-------------------------|--|
| <i>inUserOperations</i> | determine if APS is in User Operations mode<br>(boolean) |
|-------------------------|--|

---

#### **inUserOperations**

determine if APS is in User Operations mode (boolean)

Use this property to configure ophyd Devices for direct or simulated hardware. See issue #49 ([https://github.com/BCDA-APS/APS\\_BlueSky\\_tools/issues/49](https://github.com/BCDA-APS/APS_BlueSky_tools/issues/49)) for details.

EXAMPLE:

```
APS = APS_BlueSky_tools.devices.ApsMachineParametersDevice(name="APS")  
  
if APS.inUserOperations:  
    suspend_APS_current = bluesky.suspenders.SuspendFloor(APS.current, 2,  
    ↳resume_thresh=10)  
    RE.install_suspender(suspend_APS_current)  
else:  
    # use pseudo shutter controls and no current suspenders  
    pass
```

**class** APS\_BlueSky\_tools.devices.**ApsOperatorMessagesDevice**(\*args, \*\*kwargs)  
general messages from the APS main control room

**class** APS\_BlueSky\_tools.devices.**ApsPssShutter**(\*args, \*\*kwargs)  
APS PSS shutter

- APS PSS shutters have separate bit PVs for open and close

- set either bit, the shutter moves, and the bit resets a short time later
- no indication that the shutter has actually moved from the bits (see [ApsPssShutterWithStatus\(\)](#) for alternative)

EXAMPLE:

```
shutter_a = ApsPssShutter("2bma:A_shutter", name="shutter")

shutter_a.open()
shutter_a.close()

shutter_a.set("open")
shutter_a.set("close")
```

When using the shutter in a plan, be sure to use `yield from`, such as:

```
def in_a_plan(shutter):
    yield from abs_set(shutter, "open", wait=True)
    # do something
    yield from abs_set(shutter, "close", wait=True)

RE(in_a_plan(shutter_a))
```

The strings accepted by `set()` are defined in two lists: `valid_open_values` and `valid_close_values`. These lists are treated (internally to `set()`) as lower case strings.

Example, add “o” & “x” as aliases for “open” & “close”:

```
shutter_a.valid_open_values.append("o")      shutter_a.valid_close_values.append("x")      shut-
ter_a.set("o") shutter_a.set("x")

close()
request shutter to close, interactive use

open()
request shutter to open, interactive use

set(value, **kwargs)
request the shutter to open or close, BlueSky plan use

class APS_BlueSky_tools.devices.ApsPssShutterWithStatus(prefix, state_pv, *args,
**kwargs)
```

APS PSS shutter with separate status PV

- APS PSS shutters have separate bit PVs for open and close
- set either bit, the shutter moves, and the bit resets a short time later
- a separate status PV tells if the shutter is open or closed (see [ApsPssShutter\(\)](#) for alternative)

EXAMPLE:

```
A_shutter = ApsPssShutterWithStatus(
    "2bma:A_shutter",
    "PA:02BM:STA_A_FES_OPEN_PL",
    name="A_shutter")
B_shutter = ApsPssShutterWithStatus(
    "2bma:B_shutter",
    "PA:02BM:STA_B_SBS_OPEN_PL",
    name="B_shutter")
```

(continues on next page)

(continued from previous page)

```
A_shutter.open()
A_shutter.close()

or

%mov A_shutter "open"
%mov A_shutter "close"

or

A_shutter.set("open")           # MUST be "open", not "Open"
A_shutter.set("close")
```

When using the shutter in a plan, be sure to use *yield from*.

```
def in_a_plan(shutter): yield from abs_set(shutter, "open", wait=True) # do something
                                yield from abs_set(shutter, "close", wait=True)

RE(in_a_plan(A_shutter))
```

The strings accepted by *set()* are defined in attributes (*open\_str* and *close\_str*).

**close** (*timeout*=10)

**isClosed**

**isOpen**

**open** (*timeout*=10)

```
class APS_BlueSky_tools.devices.ApsUndulator(*args, **kwargs)
APS Undulator
```

EXAMPLE:

```
undulator = ApsUndulator("ID09ds:", name="undulator")
```

```
class APS_BlueSky_tools.devices.ApsUndulatorDual(*args, **kwargs)
APS Undulator with upstream and downstream controls
```

EXAMPLE:

```
undulator = ApsUndulatorDual("ID09", name="undulator")
```

note:: the trailing : in the PV prefix should be omitted

```
exception APS_BlueSky_tools.devices.AxisTunerException
Exception during execution of AxisTunerBase subclass
```

```
class APS_BlueSky_tools.devices.AxisTunerMixin(*args, **kwargs)
Mixin class to provide tuning capabilities for an axis
```

See the *TuneAxis()* example in this jupyter notebook: [https://github.com/BCDA-APS/APS\\_BlueSky\\_tools/blob/master/docs/source/resources/demo\\_tuneaxis.ipynb](https://github.com/BCDA-APS/APS_BlueSky_tools/blob/master/docs/source/resources/demo_tuneaxis.ipynb)

## HOOK METHODS

There are two hook methods (*pre\_tune\_method()*, and *post\_tune\_method()*) for callers to add additional plan parts, such as opening or closing shutters, setting detector parameters, or other actions.

Each hook method must accept a single argument: an axis object such as *EpicsMotor* or *SynAxis*, such as:

```

def my_pre_tune_hook(axis):
    yield from bps.mv(shutter, "open")
def my_post_tune_hook(axis):
    yield from bps.mv(shutter, "close")

class TunableSynAxis(AxisTunerMixin, SynAxis): pass

myaxis = TunableSynAxis(name="myaxis")
mydet = SynGauss('mydet', myaxis, 'myaxis', center=0.21, Imax=0.98e5, sigma=0.127)
myaxis.tuner = TuneAxis([mydet], myaxis)
myaxis.pre_tune_method = my_pre_tune_hook
myaxis.post_tune_method = my_post_tune_hook

RE(myaxis.tune())

```

**class** APS\_BlueSky\_tools.devices.DeviceMixinBase(\*args, \*\*kwargs)  
Base class for APS\_Bluessky\_tools Device mixin classes

**class** APS\_BlueSky\_tools.devices.DualPf4FilterBox(\*args, \*\*kwargs)  
Dual Xia PF4 filter boxes using support from synApps (using Al, Ti foils)

EXAMPLE:

```

pf4 = DualPf4FilterBox("2bmb:pf4:", name="pf4")
pf4_AlTi = DualPf4FilterBox("9idcRIO:pf4:", name="pf4_AlTi")

```

**class** APS\_BlueSky\_tools.devices.EpicsDescriptionMixin(\*args, \*\*kwargs)  
add a record's description field to a Device, such as EpicsMotor

EXAMPLE:

```

from ophyd import EpicsMotor
from APS_BlueSky_tools.devices import EpicsDescriptionMixin

class myEpicsMotor(EpicsDescriptionMixin, EpicsMotor): pass
m1 = myEpicsMotor('xxx:m1', name='m1')
print(m1.desc.value)

```

**class** APS\_BlueSky\_tools.devices.EpicsMotorDialMixin(\*args, \*\*kwargs)  
add motor record's dial coordinate fields to Device

EXAMPLE:

```

from ophyd import EpicsMotor
from APS_BlueSky_tools.devices import EpicsMotorDialMixin

class myEpicsMotor(EpicsMotorDialMixin, EpicsMotor): pass
m1 = myEpicsMotor('xxx:m1', name='m1')
print(m1.dial.read())

```

**class** APS\_BlueSky\_tools.devices.EpicsMotorLimitsMixin(\*args, \*\*kwargs)  
add motor record HLM & LLM fields & compatibility get\_lim() and set\_lim()

EXAMPLE:

```

from ophyd import EpicsMotor
from APS_BlueSky_tools.devices import EpicsMotorLimitsMixin

```

(continues on next page)

(continued from previous page)

```
class myEpicsMotor(EpicsMotorLimitsMixin, EpicsMotor): pass
m1 = myEpicsMotor('xxx:m1', name='m1')
lo = m1.get_lim(-1)
hi = m1.get_lim(1)
m1.set_lim(-25, -5)
print(m1.get_lim(-1), m1.get_lim(1))
m1.set_lim(lo, hi)
```

**get\_lim(*flag*)**

Returns the user limit of motor

- flag > 0: returns high limit
- flag < 0: returns low limit
- flag == 0: returns None

Similar with SPEC command

**set\_lim(*low, high*)**

Sets the low and high limits of motor

- No action taken if motor is moving.
- Low limit is set to lesser of (low, high)
- High limit is set to greater of (low, high)

Similar with SPEC command

```
class APS_BlueSky_tools.devices.EpicsMotorRawMixin(*args, **kwargs)
add motor record's raw coordinate fields to Device
```

EXAMPLE:

```
from ophyd import EpicsMotor
from APS_BlueSky_tools.devices import EpicsMotorRawMixin

class myEpicsMotor(EpicsMotorRawMixin, EpicsMotor): pass
m1 = myEpicsMotor('xxx:m1', name='m1')
print(m1.raw.read())
```

```
class APS_BlueSky_tools.devices.EpicsMotorServoMixin(*args, **kwargs)
add motor record's servo loop controls to Device
```

EXAMPLE:

```
from ophyd import EpicsMotor
from APS_BlueSky_tools.devices import EpicsMotorServoMixin

class myEpicsMotor(EpicsMotorServoMixin, EpicsMotor): pass
m1 = myEpicsMotor('xxx:m1', name='m1')
print(m1.servo.read())
```

```
class APS_BlueSky_tools.devices.EpicsMotorShutter(*args, **kwargs)
a shutter, implemented with an EPICS motor moved between two positions
```

EXAMPLE:

```
tomo_shutter = EpicsMotorShutter("2bma:m23", name="tomo_shutter")
tomo_shutter.closed_position = 1.0      # default
tomo_shutter.open_position = 0.0       # default
tomo_shutter.open()
tomo_shutter.close()

# or, when used in a plan
def planA():
    yield from abs_set(tomo_shutter, "open", group="O")
    yield from wait("O")
    yield from abs_set(tomo_shutter, "close", group="X")
    yield from wait("X")
def planA():
    yield from abs_set(tomo_shutter, "open", wait=True)
    yield from abs_set(tomo_shutter, "close", wait=True)
def planA():
    yield from mv(tomo_shutter, "open")
    yield from mv(tomo_shutter, "close")
```

**close()**

move motor to BEAM BLOCKED position, interactive use

**isClosed****isOpen****open()**

move motor to BEAM NOT BLOCKED position, interactive use

**set (value, \*, timeout=None, settle\_time=None)**

*set()* is like *put()*, but used in BlueSky plans

**PARAMETERS**

value : “open” or “close”

**timeout** [float, optional] Maximum time to wait. Note that *set\_and\_wait* does not support an infinite timeout.

**settle\_time: float, optional** Delay after the *set()* has completed to indicate completion to the caller

**RETURNS**

status : DeviceStatus

**class APS\_BlueSky\_tools.devices.EpicsOnOffShutter(\*args, \*\*kwargs)**

a shutter, implemented with an EPICS PV moved between two positions

Use for a shutter controlled by a single PV which takes a value for the close command and a different value for the open command. The current position is determined by comparing the value of the control with the expected open and close values.

**EXAMPLE:**

```
bit_shutter = EpicsOnOffShutter("2bma:bit1", name="bit_shutter")
bit_shutter.closed_position = 0      # default
bit_shutter.open_position = 1       # default
bit_shutter.open()
bit_shutter.close()

# or, when used in a plan
```

(continues on next page)

(continued from previous page)

```
def planA():
    yield from mv(bit_shutter, "open")
    yield from mv(bit_shutter, "close")
```

**close()**

move control to BEAM BLOCKED position, interactive use

**isClosed****isOpen****open()**

move control to BEAM NOT BLOCKED position, interactive use

**set(value, \*, timeout=None, settle\_time=None)**

*set()* is like *put()*, but used in BlueSky plans

**PARAMETERS**

value : “open” or “close”

**timeout** [float, optional] Maximum time to wait. Note that *set\_and\_wait* does not support an infinite timeout.

**settle\_time: float, optional** Delay after the *set()* has completed to indicate completion to the caller

**RETURNS**

status : DeviceStatus

**class APS\_BlueSky\_tools.devices.ProcedureRegistry(\*args, \*\*kwargs)**

Procedure Registry: run a blocking function in a thread

With many instruments, such as USAXS, there are several operating modes to be used, each with its own setup code. This ophyd Device should coordinate those modes so that the setup procedures can be called either as part of a Bluesky plan or from the command line directly. Assumes that users will write functions to setup a particular operation or operating mode. The user-written functions may not be appropriate to use in a plan directly since they might make blocking calls. The ProcedureRegistry will call the function in a thread (which is allowed to make blocking calls) and wait for the thread to complete.

It is assumed that each user-written function will not return until it is complete. .. autosummary:

```
~dir
~add
~remove
~set
~put
```

**EXAMPLE:**

Given these function definitions:

```
def clearScalerNames():
    for ch in scaler.channels.configuration_attrs:
        if ch.find(".") < 0:
            chan = scaler.channels.__getattribute__(ch)
            chan.chname.put("")

def setMyScalerNames():
    scaler.channels.chan01.chname.put("clock")
```

(continues on next page)

(continued from previous page)

```
scaler.channels.chan02.chname.put("IO")
scaler.channels.chan03.chname.put("detector")
```

create a registry and add the two functions (default name is the function name):

```
use_mode = ProcedureRegistry(name="ProcedureRegistry") use_mode.add(clearScalerNames)
use_mode.add(setMyScalerNames)
```

and then use this registry in a plan, such as this:

```
def myPlan():
    yield from bps.mv(use_mode, "setMyScalerNames")
    yield from bps.sleep(5)
    yield from bps.mv(use_mode, "clearScalerNames")
```

**add**(procedure, proc\_name=None)  
add procedure to registry

**dir**  
tuple of procedure names

**put**(value)  
replaces ophyd Device default put() behavior

**remove**(procedure)  
remove procedure from registry

**set**(proc\_name)  
run procedure in a thread, return once it is complete  
proc\_name (str) : name of registered procedure to be run

```
class APS_BlueSky_tools.devices.SimulatedApsPssShutterWithStatus(*args,
                                                               **kwargs)
```

Simulated APS PSS shutter

EXAMPLE:

```
sim = SimulatedApsPssShutterWithStatus(name="sim")
```

**close**(timeout=10)  
request the shutter to close

**get\_response\_time**()  
simulated response time for PSS status

**isClosed**  
is the shutter closed?

**isOpen**  
is the shutter open?

**open**(timeout=10)  
request the shutter to open

**set**(value, \*\*kwargs)  
set the shutter to “close” or “open”

```
APS_BlueSky_tools.devices.use_EPICS_scaler_channels(scaler)
```

configure scaler for only the channels with names assigned in EPICS

## 1.6 File Writers

BlueSky callback that writes SPEC data files

---

|   |   |
|---|---|
| <code>SpecWriterCallback([filename, auto_write])</code> | collect data from BlueSky RunEngine documents to write as SPEC data |
|---|---|

---

EXAMPLE : the `specfile_example()` writes one or more scans to a SPEC data file using a jupyter notebook.

EXAMPLE : use as BlueSky callback:

```
from APS_BlueSky_tools.filewriters import SpecWriterCallback
specwriter = SpecWriterCallback()
RE.subscribe(specwriter.receiver)
```

EXAMPLE : use as writer from Databroker:

```
from APS_BlueSky_tools.filewriters import SpecWriterCallback
specwriter = SpecWriterCallback()
for key, doc in db.get_documents(db[-1]):
    specwriter.receiver(key, doc)
print("Look at SPEC data file: "+specwriter.spec_filename)
```

EXAMPLE : use as writer from Databroker with customizations:

```
from APS_BlueSky_tools.filewriters import SpecWriterCallback

# write into file: /tmp/cerium.spec
specwriter = SpecWriterCallback(filename="/tmp/cerium.spec")
for key, doc in db.get_documents(db[-1]):
    specwriter.receiver(key, doc)

# write into file: /tmp/barium.dat
specwriter.newfile("/tmp/barium.dat")
for key, doc in db.get_documents(db["b46b63d4"]):
    specwriter.receiver(key, doc)
```

**class** `APS_BlueSky_tools.filewriters.SpecWriterCallback(filename=None, auto_write=True)`

collect data from BlueSky RunEngine documents to write as SPEC data

This gathers data from all documents and appends scan to the file when the `stop` document is received.

Parameters

**filename** [string, optional] Local, relative or absolute name of SPEC data file to be used. If `filename=None`, defaults to format of YYYYmmdd-HHMMSS.dat derived from the current system time.

**auto\_write** [boolean, optional] If True (default), `write_scan()` is called when `stop` document is received. If False, the caller is responsible for calling `write_scan()` before the next `start` document is received.

User Interface methods

---

|   |  |
|---|--|
| <code>receiver(key, document)</code>                | BlueSky callback: receive all documents for handling |
| <code>newfile([filename, reset_scan_id, RE])</code> | prepare to use a new SPEC data file                  |
| <code>usefile(filename)</code>                      | read from existing SPEC data file                    |

---

Continued on next page

Table 14 – continued from previous page

|                                      |  |
|--------------------------------------|--|
| <code>make_default_filename()</code> | generate a file name to be used as default         |
| <code>clear()</code>                 | reset all scan data defaults                       |
| <code>prepare_scan_contents()</code> | format the scan for a SPEC data file               |
| <code>write_scan()</code>            | write the most recent (completed) scan to the file |

## Internal methods

|                               |  |
|-------------------------------|--|
| <code>write_header()</code>   | write the header section of a SPEC data file |
| <code>start(doc)</code>       | handle <i>start</i> documents                |
| <code>descriptor(doc)</code>  | handle <i>descriptor</i> documents           |
| <code>event(doc)</code>       | handle <i>event</i> documents                |
| <code>bulk_events(doc)</code> | handle <i>bulk_events</i> documents          |
| <code>datum(doc)</code>       | handle <i>datum</i> documents                |
| <code>resource(doc)</code>    | handle <i>resource</i> documents             |
| <code>stop(doc)</code>        | handle <i>stop</i> documents                 |

**`bulk_events(doc)`**  
handle *bulk\_events* documents

**`clear()`**  
reset all scan data defaults

**`datum(doc)`**  
handle *datum* documents

**`descriptor(doc)`**  
handle *descriptor* documents

prepare for primary scan data, ignore any other data stream

**`event(doc)`**  
handle *event* documents

**`make_default_filename()`**  
generate a file name to be used as default

**`newfile(filename=None, reset_scan_id=False, RE=None)`**  
prepare to use a new SPEC data file  
but don't create it until we have data

**`prepare_scan_contents()`**  
format the scan for a SPEC data file

**Returns** [str] a list of lines to append to the data file

**`receiver(key, document)`**  
BlueSky callback: receive all documents for handling

**`resource(doc)`**  
handle *resource* documents

**`start(doc)`**  
handle *start* documents

**`stop(doc)`**  
handle *stop* documents

**usefile (filename)**  
 read from existing SPEC data file

**write\_header ()**  
 write the header section of a SPEC data file

**write\_scan ()**  
 write the most recent (completed) scan to the file

- creates file if not existing
- writes header if needed
- appends scan data

note: does nothing if there are no lines to be written

Example output from SpecWriterCallback ():

```

1 #F test_specdata.txt
2 #E 1510948301
3 #D Fri Nov 17 13:51:41 2017
4 #C BlueSky user = mintadmin host = mint-vm
5
6 #S 233 scan(detectors=['synthetic_pseudovoigt'], num=20, motor=['m1'], start=-1.65, ↴
7 →stop=-1.25, per_step=None)
8 #D Fri Nov 17 11:58:56 2017
9 #C Fri Nov 17 11:58:56 2017. plan_type = generator
10 #C Fri Nov 17 11:58:56 2017. uid = ddb81ac5-f3ee-4219-b047-c1196d08a5c1
11 #MD beamline_id = developer__YOUR_BEAMLINE_HERE
12 #MD login_id = mintadmin@mint-vm
13 #MD motors = ['m1']
14 #MD num_intervals = 19
15 #MD num_points = 20
16 #MD pid = 7133
17 #MD plan_pattern = linspace
18 #MD plan_pattern_args = {'start': -1.65, 'stop': -1.25, 'num': 20}
19 #MD plan_pattern_module = numpy
20 #MD proposal_id = None
21 #N 20
22 #L m1 m1_user_setpoint Epoch_float Epoch synthetic_pseudovoigt
23 -1.6500000000000001 -1.65 8.27465009689331 8 2155.6249784809206
24 -1.6288 -1.6289473684210525 8.46523666381836 8 2629.5229081466964
25 -1.608 -1.6078947368421053 8.665581226348877 9 3277.4074328018964
26 -1.5868 -1.5868421052631578 8.865738153457642 9 4246.145049452576
27 -1.5656 -1.5657894736842104 9.066259145736694 9 5825.186516381953
28 -1.5448000000000002 -1.5447368421052632 9.266754627227783 9 8803.414029867528
29 -1.5236 -1.5236842105263158 9.467074871063232 9 15501.419687691103
30 -1.5028000000000001 -1.5026315789473683 9.667330741882324 10 29570.38936784884
31 -1.4816 -1.4815789473684209 9.867793798446655 10 55562.3437459487
32 -1.4604000000000001 -1.4605263157894737 10.067811012268066 10 89519.64275090238
33 -1.4396 -1.4394736842105262 10.268356084823608 10 97008.97190269837
34 -1.4184 -1.418421052631579 10.470621824264526 10 65917.29757650592
35 -1.3972 -1.3973684210526316 10.669955730438232 11 36203.46726798266
36 -1.3764 -1.3763157894736842 10.870310306549072 11 18897.64061096024
37 -1.3552 -1.3552631578947367 11.070487976074219 11 10316.223844200193
38 -1.3344 -1.3342105263157895 11.271018743515015 11 6540.179615556269
39 -1.3132000000000001 -1.313157894736842 11.4724280834198 11 4643.555421314616
40 -1.292 -1.2921052631578946 11.673305034637451 12 3533.8582404216445
-1.2712 -1.2710526315789474 11.874176025390625 12 2809.1872596809008

```

(continues on next page)

(continued from previous page)

```

41 -1.25 -1.25 12.074703216552734 12 2285.9226305883626
42 #C Fri Nov 17 11:59:08 2017. num_events_primary = 20
43 #C Fri Nov 17 11:59:08 2017. time = 2017-11-17 11:59:08.324011
44 #C Fri Nov 17 11:59:08 2017. exit_status = success

```

## 1.7 Plans

Plans that might be useful at the APS when using BlueSky

|   |  |
|---|--|
| <code>nscan(detectors, *motor_sets[, num, ...])</code>  | Scan over n variables moved together, each in equally spaced steps.                    |
| <code>ProcedureRegistry(*args, **kwargs)</code>         | Procedure Registry   |
| <code>run_blocker_in_plan(blocker, *args[, ...])</code> | plan: run blocking function <code>blocker_(*args, **kwargs)</code> from a Bluesky plan |
| <code>run_in_thread(func)</code>                        | (decorator) run <code>func</code> in thread  |
| <code>snapshot(obj_list[, stream, md])</code>           | bluesky plan: record current values of list of ophyd signals                           |
| <code>TuneAxis(signals, axis[, signal_name])</code>     | tune an axis with a signal   |
| <code>tune_axes(axes)</code>                            | BlueSky plan to tune a list of axes in sequence  |

```

class APS_BlueSky_tools.plans.ProcedureRegistry(*args, **kwargs)
    Procedure Registry

```

**Caution:** This Device may be relocated or removed entirely in future releases. Its use is complicated and could lead to instability.

With many instruments, such as USAXS, there are several operating modes to be used, each with its own setup code. This ophyd Device should coordinate those modes so that the setup procedures can be called either as part of a Bluesky plan or from the command line directly.

Assumes that users will write functions to setup a particular operation or operating mode. The user-written functions may not be appropriate to use in a plan directly since they might make blocking calls. The ProcedureRegistry will call the function in a thread (which is allowed to make blocking calls) and wait for the thread to complete.

It is assumed that each user-written function will not return until it is complete.

|  |   |
|--|---|
| <code>dir</code>                         | tuple of procedure names                              |
| <code>add(procedure[, proc_name])</code> | add procedure to registry                             |
| <code>remove(procedure)</code>           | remove procedure from registry                        |
| <code>set(proc_name)</code>              | run procedure in a thread, return once it is complete |
| <code>put(value)</code>                  | replaces ophyd Device default put() behavior          |

EXAMPLE:

```

use_mode = ProcedureRegistry(name="use_mode")

def clearScalerNames():
    for ch in scaler.channels.configuration_attrs:

```

(continues on next page)

(continued from previous page)

```

if ch.find(".") < 0:
    chan = scaler.channels.__getattribute__(ch)
    chan.chname.put("")

def setMyScalerNames():
    scaler.channels.chan01.chname.put("clock")
    scaler.channels.chan02.chname.put("IO")
    scaler.channels.chan03.chname.put("detector")

def useMyScalerNames(): # Bluesky plan
    yield from bps.mv(
        m1, 5,
        use_mode, "clear",
    )
    yield from bps.mv(
        m1, 0,
        use_mode, "set",
    )

def demo():
    print(1)
    m1.move(5)
    print(2)
    time.sleep(2)
    print(3)
    m1.move(0)
    print(4)

use_mode.add(demo)
use_mode.add(clearScalerNames, "clear")
use_mode.add(setMyScalerNames, "set")
# use_mode.set("demo")
# use_mode.set("clear")
# RE(useMyScalerNames())

```

**add**(procedure, proc\_name=None)  
add procedure to registry

**dir**  
tuple of procedure names

**put**(value)  
replaces ophyd Device default put() behavior

**remove**(procedure)  
remove procedure from registry

**set**(proc\_name)  
run procedure in a thread, return once it is complete  
proc\_name (str) : name of registered procedure to be run

**class** APS\_BlueSky\_tools.plans.TuneAxis(signals, axis, signal\_name=None)  
tune an axis with a signal

This class provides a tuning object so that a Device or other entity may gain its own tuning process, keeping track of the particulars needed to tune this device again. For example, one could add a tuner to a motor stage:

```
motor = EpicsMotor("xxx:motor", "motor")
motor.tuner = TuneAxis([det], motor)
```

Then the motor could be tuned individually:

```
RE(motor.tuner.tune(md={"activity": "tuning"}))
```

or the `tune()` could be part of a plan with other steps.

Example:

```
tuner = TuneAxis([det], axis)
live_table = LiveTable(["axis", "det"])
RE(tuner.multi_pass_tune(width=2, num=9), live_table)
RE(tuner.tune(width=0.05, num=9), live_table)
```

Also see the jupyter notebook referenced here: [Example: TuneAxis\(\)](#).

|  |   |
|--|---|
| <code>tune([width, num, md])</code>                          | BlueSky plan to execute one pass through the current scan range |
| <code>multi_pass_tune([width, step_factor, num, ...])</code> | BlueSky plan for tuning this axis with this signal              |
| <code>peak_detected()</code>                                 | returns True if a peak was detected, otherwise False            |

**`multi_pass_tune(width=None, step_factor=None, num=None, pass_max=None, snake=None, md=None)`**

BlueSky plan for tuning this axis with this signal

Execute multiple passes to refine the centroid determination. Each subsequent pass will reduce the width of scan by `step_factor`. If `snake=True` then the scan direction will reverse with each subsequent pass.

#### PARAMETERS

**`width`** [float] width of the tuning scan in the units of `self.axis` Default value in `self.width` (initially 1)

**`num`** [int] number of steps Default value in `self.num` (initially 10)

**`step_factor`** [float] This reduces the width of the next tuning scan by the given factor. Default value in `self.step_factor` (initially 4)

**`pass_max`** [int] Maximum number of passes to be executed (avoids runaway scans when a centroid is not found). Default value in `self.pass_max` (initially 10)

**`snake`** [bool] If `True`, reverse scan direction on next pass. Default value in `self.snake` (initially `True`)

**`md`** [dict, optional] metadata

**`peak_detected()`**

returns True if a peak was detected, otherwise False

The default algorithm identifies a peak when the maximum value is four times the minimum value. Change this routine by subclassing `TuneAxis` and override `peak_detected()`.

**`tune(width=None, num=None, md=None)`**

BlueSky plan to execute one pass through the current scan range

Scan `self.axis` centered about current position from `-width/2` to `+width/2` with `num` observations. If a peak was detected (default check is that `max >= 4*min`), then set `self.tune_ok = True`.

## PARAMETERS

**width** [float] width of the tuning scan in the units of `self.axis`. Default value in `self.width` (initially 1)

**num** [int] number of steps Default value in `self.num` (initially 10)

**md** [dict, optional] metadata

`APS_BlueSky_tools.plans.nscan(detectors, *motor_sets, num=11, per_step=None, md=None)`

Scan over n variables moved together, each in equally spaced steps.

## PARAMETERS

**detectors** [list] list of ‘readable’ objects

**motor\_sets** [list] sequence of one or more groups of: motor, start, finish

**motor** [object] any ‘settable’ object (motor, temp controller, etc.)

**start** [float] starting position of motor

**finish** [float] ending position of motor

**num** [int] number of steps (default = 11)

**per\_step** [callable, optional] hook for customizing action of inner loop (messages per step) Expected signature:  
`f(detectors, step_cache, pos_cache)`

**md** [dict, optional] metadata

See the `nscan()` example in a Jupyter notebook: [https://github.com/BCDA-APS/APS\\_BlueSky\\_tools/blob/master/docs/source/resources/demo\\_nscan.ipynb](https://github.com/BCDA-APS/APS_BlueSky_tools/blob/master/docs/source/resources/demo_nscan.ipynb)

`APS_BlueSky_tools.plans.run_blocker_in_plan(blocker, *args, _poll_s_=0.01, _timeout_s_=None, **kwargs)`

plan: run blocking function `blocker_(*args, **kwargs)` from a Bluesky plan

## PARAMETERS

**blocker** [func] function object to be called in a Bluesky plan

**\_poll\_s\_** [float] sleep interval in loop while waiting for completion (default: 0.01)

**\_timeout\_s\_** [float] maximum time for completion (default: *None* which means no timeout)

Example: use `time.sleep` as blocking function:

```
RE(run_blocker_in_plan(time.sleep, 2.14))
```

Example: in a plan, use `time.sleep` as blocking function:

```
def my_sleep(t=1.0):
    yield from run_blocker_in_plan(time.sleep, t)

RE(my_sleep())
```

`APS_BlueSky_tools.plans.run_in_thread(func)`  
 (decorator) run `func` in thread

USAGE:

```
@run_in_thread
def progress_reporting():
    logger.debug("progress_reporting is starting")
```

(continues on next page)

(continued from previous page)

```
# ...  
  
#...  
progress_reporting()    # runs in separate thread  
#...
```

APS\_BlueSky\_tools.plans.**snapshot** (*obj\_list*, *stream*=’primary’, *md*=None)  
bluesky plan: record current values of list of ophyd signals

#### PARAMETERS

**obj\_list** [list] list of ophyd Signal or EpicsSignal objects

**stream** [str] document stream, default: “primary”

**md** [dict] metadata

APS\_BlueSky\_tools.plans.**tune\_axes** (*axes*)  
BlueSky plan to tune a list of axes in sequence

#### EXAMPLE

Sequentially, tune a list of preconfigured axes:

```
RE(tune_axes([mr, m2r, ar, a2r]))
```

## 1.8 Signals

(ophyd) Signals that might be useful at the APS using Bluesky

---

*SynPseudoVoigt*(*name*, *motor*, *motor\_field*[, …])      Evaluate a point on a pseudo-Voigt based on the value of a motor.

---

**class** APS\_BlueSky\_tools.signals.**SynPseudoVoigt** (*name*, *motor*, *motor\_field*, *center*=0, *eta*=0.5, *scale*=1, *sigma*=1, *bkg*=0, *noise*=None, *noise\_multiplier*=1, *\*\*kwargs*)

Evaluate a point on a pseudo-Voigt based on the value of a motor.

Provides a signal to be measured. Acts like a detector.

See [https://en.wikipedia.org/wiki/Voigt\\_profile](https://en.wikipedia.org/wiki/Voigt_profile)

#### PARAMETERS

**name** [str] name of detector signal

**motor** [*Mover*] The independent coordinate

**motor\_field** [str] name of *Mover* field

**center** [float, optional] location of maximum value, default=0

**eta** [float, optional]  $0 \leq \text{eta} < 1.0$ : Lorentzian fraction, default=0.5

**scale** [float, optional] scale  $\geq 1$  : scale factor, default=1

**sigma** [float, optional] sigma  $> 0$  : width, default=1

**bkg** [float, optional] bkg  $\geq 0$  : constant background, default=0

**noise** [{‘poisson’, ‘uniform’, None}] Add noise to the result.

**noise\_multiplier** [float] Only relevant for ‘uniform’ noise. Multiply the random amount of noise by ‘noise\_multiplier’

EXAMPLE

```
from APS_BlueSky_tools.signals import SynPseudoVoigt
motor = Mover('motor', {'motor': lambda x: x}, {'x': 0})
det = SynPseudoVoigt('det', motor, 'motor',
    center=0, eta=0.5, scale=1, sigma=1, bkg=0)
```

EXAMPLE

```
import numpy as np
from APS_BlueSky_tools.signals import SynPseudoVoigt
synthetic_pseudovoigt = SynPseudoVoigt(
    'synthetic_pseudovoigt', m1, 'm1',
    center=-1.5 + 0.5*np.random.uniform(),
    eta=0.2 + 0.5*np.random.uniform(),
    sigma=0.001 + 0.05*np.random.uniform(),
    scale=1e5,
    bkg=0.01*np.random.uniform())

# RE(bp.scan([synthetic_pseudovoigt], m1, -2, 0, 219))
```

## 1.9 Suspenders

(bluesky) custom support for pausing a running plan

---

|   |                   |
|---|-------------------|
| <i>SuspendWhenChanged(signal, *[...])</i> | Bluesky suspender |
|---|-------------------|

---

```
class APS_BlueSky_tools.suspenders.SuspendWhenChanged(signal, *, expected_value=None, allow_resume=False, sleep=0, pre_plan=None, post_plan=None, tripped_message='', **kwargs)
```

Bluesky suspender

Suspend when the monitored value deviates from the expected. Only resume if allowed AND when monitored equals expected. Default expected value is current value when object is created.

USAGE:

```
# pause if this value changes in our session
# note: this suspender is designed to require Bluesky restart if value changes
suspend_instrument_in_use = SuspendWhenChanged(instrument_in_use)
RE.install_suspender(suspend_instrument_in_use)
```

## 1.10 Utilities

Various utilities

|   |   |
|---|---|
| <code>connect_pvlist(pvlist[, wait, timeout, ...])</code>     | given a list of EPICS PV names, return a dictionary of EpicsSignal objects              |
| <code>EmailNotifications([sender])</code>                     | send email notifications when requested   |
| <code>ExcelDatabaseFileBase()</code>                          | base class: read-only support for Excel files, treat them like databases                |
| <code>ExcelDatabaseFileGeneric(filename[, labels_row])</code> | Generic (read-only) handling of Excel spreadsheet-as-database                           |
| <code>ipython_profile_name()</code>                           | return the name of the current ipython profile or <i>None</i>                           |
| <code>print_snapshot_list(db, **search_criteria)</code>       | print (stdout) a list of all snapshots in the databroker                                |
| <code>text_encode(source)</code>                              | encode source using the default codepoint   |
| <code>to_unicode_or_bust(obj[, encoding])</code>              | from: <a href="http://farmdev.com/talks/unicode/">http://farmdev.com/talks/unicode/</a> |
| <code>unix_cmd(command_list)</code>                           | run a UNIX command, returns (stdout, stderr)  |

**class** APS\_BlueSky\_tools.utils.**EmailNotifications** (*sender=None*)  
send email notifications when requested

use default OS mail utility (so no credentials needed)

**send** (*subject, message*)  
send message to all addresses

**class** APS\_BlueSky\_tools.utils.**ExcelDatabaseFileBase**  
base class: read-only support for Excel files, treat them like databases

#### EXAMPLE

Show how to read an Excel file where one of the columns contains a unique key. This allows for random access to each row of data by use of the *key*.

```
class ExhibitorsDB(ExcelDatabaseFileBase):
    """
    content for Exhibitors, vendors, and Sponsors from the Excel file
    """

    EXCEL_FILE = os.path.join("resources", "exhibitors.xlsx")
    LABELS_ROW = 2

    def handle_single_entry(self, entry):
        '''any special handling for a row from the Excel file'''
        pass

    def handleExcelRowEntry(self, entry):
        '''identify the unique key for this entry (row of the Excel file)'''
        key = entry["Name"]
        self.db[key] = entry
```

**class** APS\_BlueSky\_tools.utils.**ExcelDatabaseFileGeneric** (*filename, labels\_row=3*)  
Generic (read-only) handling of Excel spreadsheet-as-database

Table labels are given on Excel row N, *self.labels\_row* = N-1

**handleExcelRowEntry** (*entry*)  
use row number as the unique key

APS\_BlueSky\_tools.utils.**connect\_pvlist** (*pvlist, wait=True, timeout=2, poll\_interval=0.1*)  
given a list of EPICS PV names, return a dictionary of EpicsSignal objects

#### PARAMETERS

**pvlist** [list(str)] list of EPICS PV names

**wait** [bool] should wait for EpicsSignal objects to connect, default: True  
**timeout** [float] maximum time to wait for PV connections, seconds, default: 2.0  
**poll\_interval** [float] time to sleep between checks for PV connections, seconds, default: 0.1

APS\_BlueSky\_tools.utils.**ipython\_profile\_name()**  
return the name of the current ipython profile or *None*

Example (add to default RunEngine metadata):

```
RE.md['ipython_profile'] = str(ipython_profile_name())
print("using profile: " + RE.md['ipython_profile'])
```

APS\_BlueSky\_tools.utils.**print\_snapshot\_list**(db, \*\*search\_criteria)  
print (stdout) a list of all snapshots in the databroker

USAGE:

```
print_snapshot_list(db, )
print_snapshot_list(db, purpose="this is an example")
print_snapshot_list(db, since="2018-12-21", until="2019")
```

EXAMPLE:

```
In [16]: from APS_BlueSky_tools.utils import print_snapshot_list
...: from APS_BlueSky_tools.callbacks import SnapshotReport
...: print_snapshot_list(db, since="2018-12-21", until="2019")
...:
=====
# uid      date/time            purpose
=====
0 d7831dae 2018-12-21 11:39:52.956904 this is an example
1 5049029d 2018-12-21 11:39:30.062463 this is an example
2 588e0149 2018-12-21 11:38:43.153055 this is an example
=====

In [17]: SnapshotReport().print_report(db["5049029d"])

=====
snapshot: 2018-12-21 11:39:30.062463
=====

example: example 2
hints: {}
iso8601: 2018-12-21 11:39:30.062463
look: can snapshot text and arrays too
note: no commas in metadata
plan_description: archive snapshot of ophyd Signals (usually EPICS PVs)
plan_name: snapshot
plan_type: generator
purpose: this is an example
scan_id: 1
software_versions: {
    'python':
        '''3.6.2 |Continuum Analytics, Inc.| (default, Jul 20 2017, 13:51:32)
[GCC 4.4.7 20120313 (Red Hat 4.4.7-1)]''',
    'PyEpics': '3.3.1',
    'bluesky': '1.4.1',
    'ophyd': '1.3.0',
```

(continues on next page)

(continued from previous page)

```
'databroker': '0.11.3',
'APS_Bluessky_Tools': '0.0.38'
}
time: 1545413970.063167
uid: 5049029d-075c-453c-96d2-55431273852b

=====
timestamp          source   name      value
=====
2018-12-20 18:24:34.220028 PV      compress    [0.1, 0.2, 0.3]
2018-12-13 14:49:53.121188 PV      gov:HOSTNAME otz.aps.anl.gov
2018-12-21 11:39:24.268148 PV      gov:IOC_CPU_LOAD 0.22522317161410768
2018-12-21 11:39:24.268151 PV      gov:SYS_CPU_LOAD 9.109026666525944
2018-12-21 11:39:30.017643 PV      gov:iso8601    2018-12-21T11:39:30
2018-12-13 14:49:53.135016 PV      otz:HOSTNAME otz.aps.anl.gov
2018-12-21 11:39:27.705304 PV      otz:IOC_CPU_LOAD 0.1251210270549924
2018-12-21 11:39:27.705301 PV      otz:SYS_CPU_LOAD 11.611234438304471
2018-12-21 11:39:30.030321 PV      otz:iso8601    2018-12-21T11:39:30
=====

exit_status: success
num_events: {'primary': 1}
run_start: 5049029d-075c-453c-96d2-55431273852b
time: 1545413970.102147
uid: 6c1b2100-1ef6-404d-943e-405da9ada882
```

`APS_BlueSky_tools.utils.text_encode(source)`  
encode source using the default codepoint

`APS_BlueSky_tools.utils.to_unicode_or_bust(obj, encoding='utf-8')`  
from: <http://farmdev.com/talks/unicode/>

`APS_BlueSky_tools.utils.unix_cmd(command_list)`  
run a UNIX command, returns (stdout, stderr)

## 1.11 synApps busy record

see the synApps busy module support: <https://github.com/epics-modules/busy>

Ophyd support for the EPICS busy record

Public Structures

---

`busyRecord(*args, **kwargs)`

---

`class APS_BlueSky_tools.synApps_ophyd.busy.busyRecord(*args, **kwargs)`

## 1.12 synApps sscan record

see the synApps sscan module support: <https://github.com/epics-modules/sscan>

Ophyd support for the EPICS synApps sscan record

EXAMPLE

```
import APS_BlueSky_tools.synApps_ophyd
scans = APS_BlueSky_tools.synApps_ophyd.sscanDevice("xxx:", name="scans")
```

## Public Structures

|   |   |
|---|---|
| <code>sscanRecord(*args, **kwargs)</code> | EPICS synApps sscan record: used as \$(P):scan(N)       |
| <code>sscanDevice(*args, **kwargs)</code> | synApps XXX IOC setup of sscan records: \$(P):scan\$(N) |

## Private Structures

|   |   |
|---|---|
| <code>sscanPositioner(prefix, num, **kwargs)</code> | positioner of an EPICS sscan record       |
| <code>sscanDetector(prefix, num, **kwargs)</code>   | detector of an EPICS sscan record         |
| <code>sscanTrigger(prefix, num, **kwargs)</code>    | detector trigger of an EPICS sscan record |

```
class APS_BlueSky_tools.synApps_ophyd.sscan.sscanRecord(*args, **kwargs)
EPICS synApps sscan record: used as $(P):scan(N)

reset()
    set all fields to default values

set(value, **kwargs)
    interface to use bps.mv()

class APS_BlueSky_tools.synApps_ophyd.sscan.sscanDevice(*args, **kwargs)
synApps XXX IOC setup of sscan records: $(P):scan$(N)

reset()
    set all fields to default values
```

## 1.13 synApps swait record

The swait record is part of the calc module: <https://htmlpreview.github.io/?https://raw.githubusercontent.com/epics-modules/calc/R3-6-1/documentation/swaitRecord.html>

see the synApps calc module support: <https://github.com/epics-modules/calc>

Ophyd support for the EPICS synApps swait record

EXAMPLES:;

```
import APS_BlueSky_tools.synApps_ophyd
calcs = APS_BlueSky_tools.synApps_ophyd.userCalcsDevice("xxx:", name="calcs")

calc1 = calcs.calc1
APS_BlueSky_tools.synApps_ophyd.swait_setup_random_number(calc1)

APS_BlueSky_tools.synApps_ophyd.swait_setup_incrementeer(calcs.calc2)

calc1.reset()
```

|  |   |
|--|---|
| <code>swaitRecord(*args, **kwargs)</code>                      | synApps swait record: used as \$(P):userCalc\$(N)       |
| <code>userCalcsDevice(*args, **kwargs)</code>                  | synApps XXX IOC setup of userCalcs: \$(P):userCalc\$(N) |
| <code>swait_setup_random_number(swait, **kw)</code>            | setup swait record to generate random numbers           |
| <code>swait_setup_gaussian(swait, motor[, center, ...])</code> | setup swait for noisy Gaussian                          |

Continued on next page

Table 25 – continued from previous page

|   |   |
|---|---|
| <code>swait_setup_lorentzian(swait, motor[, ...])</code>  | setup swait record for noisy Lorentzian |
| <code>swait_setup_incremter(swait[, scan, limit])</code>  | setup swait record as an incremter      |
| <hr/>   |   |
| <code>class APS_BlueSky_tools.synApps_ophyd.swait.swaitRecord(*args, **kwargs)</code>   |   |
| synApps swait record: used as \$(P):userCalc\$(N)   |   |
| <code>reset()</code>  |   |
| set all fields to default values  |   |
| <code>class APS_BlueSky_tools.synApps_ophyd.swait.userCalcsDevice(*args, **kwargs)</code>                                       |   |
| synApps XXX IOC setup of userCalcs: \$(P):userCalc\$(N)   |   |
| <code>reset()</code>  |   |
| set all fields to default values  |   |
| <code>APS_BlueSky_tools.synApps_ophyd.swait.swait_setup_random_number(swait, **kw)</code>                                       |   |
| setup swait record to generate random numbers   |   |
| <code>APS_BlueSky_tools.synApps_ophyd.swait.swait_setup_gaussian(swait, motor, center=0, width=1, scale=1, noise=0.05)</code>   |   |
| setup swait for noisy Gaussian  |   |
| <code>APS_BlueSky_tools.synApps_ophyd.swait.swait_setup_lorentzian(swait, motor, center=0, width=1, scale=1, noise=0.05)</code> |   |
| setup swait record for noisy Lorentzian   |   |
| <code>APS_BlueSky_tools.synApps_ophyd.swait.swait_setup_incremter(swait, scan=None, limit=100000)</code>                        |   |
| setup swait record as an incremter  |   |

## CHAPTER 2

---

### Indices and tables

---

- genindex
- modindex
- search



---

## Python Module Index

---

### a

APS\_BlueSky\_tools.callbacks, 9  
APS\_BlueSky\_tools.devices, 10  
APS\_BlueSky\_tools.examples, 9  
APS\_BlueSky\_tools.filewriters, 23  
APS\_BlueSky\_tools.plans, 26  
APS\_BlueSky\_tools.signals, 30  
APS\_BlueSky\_tools.snapshot, 6  
APS\_BlueSky\_tools.suspenders, 31  
APS\_BlueSky\_tools.synApps\_ophyd.busy,  
    34  
APS\_BlueSky\_tools.synApps\_ophyd.sscan,  
    34  
APS\_BlueSky\_tools.synApps\_ophyd.swait,  
    35  
APS\_BlueSky\_tools.utils, 31



---

## Index

---

### A

AD\_EpicsHdf5FileName (class in *APS\_BlueSky\_tools.devices*), 12  
AD\_setup\_FrameType () (in module *APS\_BlueSky\_tools.devices*), 14  
AD\_warmed\_up () (in module *APS\_BlueSky\_tools.devices*), 14  
add () (*APS\_BlueSky\_tools.devices.ProcedureRegistry* method), 22  
add () (*APS\_BlueSky\_tools.plans.ProcedureRegistry* method), 27  
*APS\_BlueSky\_tools.callbacks* (module), 9  
*APS\_BlueSky\_tools.devices* (module), 10  
*APS\_BlueSky\_tools.examples* (module), 9  
*APS\_BlueSky\_tools.filewriters* (module), 23  
*APS\_BlueSky\_tools.plans* (module), 26  
*APS\_BlueSky\_tools.signals* (module), 30  
*APS\_BlueSky\_tools.snapshot* (module), 6  
*APS\_BlueSky\_tools.suspenders* (module), 31  
*APS\_BlueSky\_tools.synApps\_ophyd.busy* (module), 34  
*APS\_BlueSky\_tools.synApps\_ophyd.sscan* (module), 34  
*APS\_BlueSky\_tools.synApps\_ophyd.swait* (module), 35  
*APS\_BlueSky\_tools.utils* (module), 31  
ApsBssUserInfoDevice (class in *APS\_BlueSky\_tools.devices*), 14  
ApsMachineParametersDevice (class in *APS\_BlueSky\_tools.devices*), 15  
ApsOperatorMessagesDevice (class in *APS\_BlueSky\_tools.devices*), 15  
ApsPssShutter (class in *APS\_BlueSky\_tools.devices*), 15  
ApsPssShutterWithStatus (class in *APS\_BlueSky\_tools.devices*), 16  
ApsUndulator (class in *APS\_BlueSky\_tools.devices*), 17  
ApsUndulatorDual (class in *APS\_BlueSky\_tools.devices*), 18

*APS\_BlueSky\_tools.devices*), 17  
in *AxisTunerException*, 17  
AxisTunerMixin (class in *APS\_BlueSky\_tools.devices*), 17

### B

bulk\_events () (*APS\_BlueSky\_tools.filewriters.SpecWriterCallback* method), 24  
busyRecord (class in *APS\_BlueSky\_tools.synApps\_ophyd.busy*), 34

### C

clear () (*APS\_BlueSky\_tools.filewriters.SpecWriterCallback* method), 24  
close () (*APS\_BlueSky\_tools.devices.ApsPssShutter* method), 16  
close () (*APS\_BlueSky\_tools.devices.ApsPssShutterWithStatus* method), 17  
close () (*APS\_BlueSky\_tools.devices.EpicsMotorShutter* method), 20  
close () (*APS\_BlueSky\_tools.devices.EpicsOnOffShutter* method), 21  
close () (*APS\_BlueSky\_tools.devices.SimulatedApsPssShutterWithStatus* method), 22  
connect\_pvlist () (in module *APS\_BlueSky\_tools.utils*), 32

### D

datum () (*APS\_BlueSky\_tools.filewriters.SpecWriterCallback* method), 24  
descriptor () (*APS\_BlueSky\_tools.callbacks.SnapshotReport* method), 10  
descriptor () (*APS\_BlueSky\_tools.filewriters.SpecWriterCallback* method), 24  
DeviceMixinBase (class in *APS\_BlueSky\_tools.devices*), 18  
dir (*APS\_BlueSky\_tools.devices.ProcedureRegistry* attribute), 22

```

dir (APS_BlueSky_tools.plans.ProcedureRegistry attribute), 27
document_contents_callback () (in module APS_BlueSky_tools.callbacks), 10
DocumentCollectorCallback (class in APS_BlueSky_tools.callbacks), 10
DualPf4FilterBox (class in APS_BlueSky_tools.devices), 18

E
EmailNotifications (class in APS_BlueSky_tools.utils), 32
EpicsDescriptionMixin (class in APS_BlueSky_tools.devices), 18
EpicsMotorDialMixin (class in APS_BlueSky_tools.devices), 18
EpicsMotorLimitsMixin (class in APS_BlueSky_tools.devices), 18
EpicsMotorRawMixin (class in APS_BlueSky_tools.devices), 19
EpicsMotorServoMixin (class in APS_BlueSky_tools.devices), 19
EpicsMotorShutter (class in APS_BlueSky_tools.devices), 19
EpicsOnOffShutter (class in APS_BlueSky_tools.devices), 20
event () (APS_BlueSky_tools.filewriters.SpecWriterCallback method), 24
ExcelDatabaseFileBase (class in APS_BlueSky_tools.utils), 32
ExcelDatabaseFileGeneric (class in APS_BlueSky_tools.utils), 32

G
generate_datum () (APS_BlueSky_tools.devices.AD_EpicsHdf5FileName module in module APS_BlueSky_tools.plans), 29
method), 14
get_args () (in module APS_BlueSky_tools.snapshot), 7
get_frames_per_point () (APS_BlueSky_tools.devices.AD_EpicsHdf5FileName method), 14
get_lim () (APS_BlueSky_tools.devices.EpicsMotorLimitsMixin method), 19
get_response_time () (APS_BlueSky_tools.devices.SimulatedApsPssShutterWithStatus method), 22

H
handleExcelRowEntry () (APS_BlueSky_tools.utils.ExcelDatabaseFileGeneric method), 32

I
inUserOperations (APS_BlueSky_tools.devices.ApsMachineParametersDevice

P
attribute), 15
ipython_profile_name () (in module APS_BlueSky_tools.utils), 33
isClosed (APS_BlueSky_tools.devices.ApsPssShutterWithStatus attribute), 17
isClosed (APS_BlueSky_tools.devices.EpicsMotorShutter attribute), 20
isClosed (APS_BlueSky_tools.devices.EpicsOnOffShutter attribute), 21
isClosed (APS_BlueSky_tools.devices.SimulatedApsPssShutterWithStatus attribute), 22
isOpen (APS_BlueSky_tools.devices.ApsPssShutterWithStatus attribute), 17
isOpen (APS_BlueSky_tools.devices.EpicsMotorShutter attribute), 20
isOpen (APS_BlueSky_tools.devices.EpicsOnOffShutter attribute), 21
isOpen (APS_BlueSky_tools.devices.SimulatedApsPssShutterWithStatus attribute), 22

M
main () (in module APS_BlueSky_tools.examples), 9
make_default_filename () (APS_BlueSky_tools.filewriters.SpecWriterCallback method), 24
make_filename () (APS_BlueSky_tools.devices.AD_EpicsHdf5FileName method), 14
multi_pass_tune () (APS_BlueSky_tools.plans.TuneAxis method), 28

N
newfile () (APS_BlueSky_tools.filewriters.SpecWriterCallback method), 24

O
open () (APS_BlueSky_tools.devices.ApsPssShutter method), 16
open () (APS_BlueSky_tools.devices.ApsPssShutterWithStatus method), 17
open () (APS_BlueSky_tools.devices.EpicsMotorShutter method), 20
open () (APS_BlueSky_tools.devices.EpicsOnOffShutter method), 21
open () (APS_BlueSky_tools.devices.SimulatedApsPssShutterWithStatus method), 22

P
peak_detected () (APS_BlueSky_tools.plans.TuneAxis method), 28
plan_catalog () (in module APS_BlueSky_tools.examples), 9

```

```

prepare_scan_contents()                               set()      (APS_BlueSky_tools.plans.ProcedureRegistry
    (APS_BlueSky_tools.filewriters.SpecWriterCallback   method), 27
    method), 24
print_report() (APS_BlueSky_tools.callbacks.SnapshotReport method), 35
    method), 10
print_snapshot_list() (in module                 set_lim() (APS_BlueSky_tools.devices.EpicsMotorLimitsMixin
    APS_BlueSky_tools.utils), 33                      method), 19
ProcedureRegistry (class in module               SimulatedApsPssShutterWithStatus (class in
    APS_BlueSky_tools.devices), 21                  APS_BlueSky_tools.devices), 22
ProcedureRegistry (class in module               snapshot() (in module APS_BlueSky_tools.plans), 30
    APS_BlueSky_tools.plans), 26
put() (APS_BlueSky_tools.devices.ProcedureRegistry SnapshotReport (class in
    method), 22                                     APS_BlueSky_tools.callbacks), 10
put() (APS_BlueSky_tools.plans.ProcedureRegistry specfile_example() (in module
    method), 27                                       APS_BlueSky_tools.examples), 9
                                                SpecWriterCallback (class in
                                                APS_BlueSky_tools.filewriters), 23
R
receiver() (APS_BlueSky_tools.callbacks.DocumentCollectorCallback receiver() (class in
    method), 10                                     APS_BlueSky_tools.synApps_ophyd.sscan),
receiver() (APS_BlueSky_tools.filewriters.SpecWriterCallback 35
    method), 24
remove() (APS_BlueSky_tools.devices.ProcedureRegistry sscanRecord (class in
    method), 22                                     APS_BlueSky_tools.synApps_ophyd.sscan),
remove() (APS_BlueSky_tools.plans.ProcedureRegistry stage() (APS_BlueSky_tools.devices.AD_EpicsHdf5FileName
    method), 27                                     method), 14
reset() (APS_BlueSky_tools.synApps_ophyd.sscan.sscanDevice) (APS_BlueSky_tools.filewriters.SpecWriterCallback
    method), 35                                     method), 24
reset() (APS_BlueSky_tools.synApps_ophyd.sscanRecord) (APS_BlueSky_tools.filewriters.SpecWriterCallback
    method), 35                                     method), 24
reset() (APS_BlueSky_tools.synApps_ophyd.swait.swaitRecord) (APS_BlueSky_tools.susppers.SuspperWhenChanged
    method), 36                                     (class in
                                                APS_BlueSky_tools.susppers), 31
reset() (APS_BlueSky_tools.synApps_ophyd.swait.userCallsDevice) (APS_BlueSky_tools.susppers.SuspperSetupGaussian
    method), 36                                     (in module
                                                APS_BlueSky_tools.synApps_ophyd.swait),
resource() (APS_BlueSky_tools.filewriters.SpecWriterCallback 36
    method), 24
run_blocker_in_plan() (in module                 swait_SetupLorentzian() (in module
    APS_BlueSky_tools.plans), 29                      APS_BlueSky_tools.synApps_ophyd.swait),
run_in_thread() (in module                     swait_SetupRandomNumber() (in module
    APS_BlueSky_tools.plans), 29                      APS_BlueSky_tools.synApps_ophyd.swait),
                                                36
S
send() (APS_BlueSky_tools.utils.EmailNotifications set()      (APS_BlueSky_tools.synApps_ophyd.swait),
    method), 32                                     method), 36
set() (APS_BlueSky_tools.devices.ApsPssShutter swaitRecord (class in
    method), 16                                     APS_BlueSky_tools.synApps_ophyd.swait),
set() (APS_BlueSky_tools.devices.EpicsMotorShutter SynPseudoVoigt (class in
    method), 20                                     APS_BlueSky_tools.signals), 30
set() (APS_BlueSky_tools.devices.EpicsOnOffShutter text_encode() (in module APS_BlueSky_tools.utils),
    method), 21                                     34
set() (APS_BlueSky_tools.devices.ProcedureRegistry +WVNS encode_or_bust() (in module
    method), 22                                     APS_BlueSky_tools.utils), 34
set() (APS_BlueSky_tools.devices.SimulatedApsPssShutterWithStatus
    method), 22

```

tune () (*APS\_BlueSky\_tools.plans.TuneAxis method*),  
    28  
tune\_axes () (*in module APS\_BlueSky\_tools.plans*),  
    30  
TuneAxis (*class in APS\_BlueSky\_tools.plans*), 27

## U

unix\_cmd () (*in module APS\_BlueSky\_tools.utils*), 34  
use\_EPICS\_scaler\_channels () (*in module  
    APS\_BlueSky\_tools.devices*), 22  
usefile () (*APS\_BlueSky\_tools.filewriters.SpecWriterCallback  
    method*), 24  
userCalcsDevice    (*class  
    in  
        APS\_BlueSky\_tools.synApps\_ophyd.swait*),  
  36

## W

write\_header () (*APS\_BlueSky\_tools.filewriters.SpecWriterCallback  
    method*), 25  
write\_scan () (*APS\_BlueSky\_tools.filewriters.SpecWriterCallback  
    method*), 25